Blood?

Is it blood?



Benzidine test



No longer used



Toxic

bre Analysis

1 Benzidine test

2 Kastle-Meyer test



Phenolphthalein

bre Analysis

1 Benzidine test

Kastle-Meyer test

3 Luminol test

ibre Analysis


Based on peroxidase properties of blood

$$H_2O_2 \rightarrow H_2O + O$$

ibre Analysis



Kastle-Meyer Test




Acid form

Colourless

DNA

Kastle-Meyer Test




Alkali form

Pink

DNA

Kastle-Meyer Test



Alkali form


Pink

Reduced form

Colourless

DNA

Kastle-Meyer Test



Kastle-Meyer Test

Suspected bloodstain

- 1 Add reduced phenolphthalein
- 2 Add H_2O_2
 - → Turn pink → Most likely blood

bre Analysis

Kastle-Meyer Test

Powerful enzymes in potatoes

- Catalyst for reaction
- Turn pink
- "False positive"